

IviumBoost10012



To increase the maximum current of your Ivium potentiostat

- Full potentiostat/galvanostat compliance
- Key capability: # maximum $\pm 100A$ @ $\pm 12V$
FRA $10\mu Hz$ - $100kHz$
- All electrochemical techniques available
- Bandwidth $> 100kHz$
- Rise time $< 50 \mu s$
- Power: $1500W$

For high power applications

- Batteries
- Fuel cells
- Solar cells
- Supercaps
- Electrolysis
- Corrosion
- Deposition



The Netherlands:

IviumTechnologies B.V.
De Zaaile 11
5612AJ Eindhoven
The Netherlands

tel. +31 40 2390600
fax. +31 40 2390601
e-mail info@ivium.eu
www.ivium.com

U.S.A.:

IviumTechnologies USA
961687 Gateway Blvd., Suite 201D
Fernandina Beach, FL 32034

tel. 800-303-3885 (toll free)
904-310-9060 (office)
fax. 904-310-9068
e-mail pete@ivium.us
www.ivium.us