




# ZView® The World's Most Popular Impedance Analysis Software



- ✓ Powerful equivalent circuit modeling
- ✓ Instant fitting with common circuits
- ✓ Data processing & graphing
- ✓ Batch file processing
- ✓ Kramers-Kronig (K-K) testing



- Modify data for area-normalization, component subtraction or addition
- Generate publication-quality graphs on a variety of axes
- Analyze and fit data using built-in or user-defined Equivalent Circuits Models
  - 9 circuit element types
  - 20 distributed element types
  - >50 parameters per model
  - >2,000 data points per spectrum

## ZView Supports:

**Scribner:** ZPlot, Z60, ZPlot-Lab, MultiStat, FuelCell, MMA, MTS4

**Solartron:** 1296, SMaRT, CellTest, ModulabCSV

**PAR:** 378, 388, 398, PowerSuite Export, Versatstat3

**Other Formats\*:** ACM Instruments, Agilent, AutoLab, BioLogic, Boukamp, CH Instruments, Gamry, Gill AC, HP, Ivium, Maccor, Newtons4th, NovoControl, PalmSens, Radiometer, Zahner, and User-Defined

\*Contact Scribner Associates for further compatibility details

